

Adjustable Pencil Hardness Tester

BGD 505 Adjustable Pencil Hardness Tester is designed according with ASTM D 3363, ISO 15184. Operator can obtain three different weights (500g, 765g and 1000g) through selecting different loading weight.



Main Technical Parameters:

- ★ Three points are touched on the tested surface (roller, pencil core)
- ★ The angle between pencil and the tested surface is 45°
- ★ Built-in spirit bubble, convenient to check level when tester is working
- ★ Pressure of pencil lead is 500g, 765g, 1000g

★ Ordering Information:

BGD 505---- Adjustable Pencil Hardness Tester