

Flash Point Tester

Flash Point Tester is designed and produced according to ISO 2719, ISO 2592, ASTM D93, D65, ASTM D 6450 etc standards. It is suitable to test flash and fire point (Except fuel oil and petroleum products whose open flash point is below 79°C)

Basic Type

◆ The design of firing adopt patent technology, keeping a stable fire through controlling the flow speed for flammable liquid. Operator can adjust the flow rate hence avoiding some dangerous and inconvenient factors which would be happened by electronic or LNG firing .

- ◆ Clear and beautiful LCD screen
- ◆ Small size, light weight and reasonable structure.



BGD 235



BGD 238

Main Technical Parameters:

- ★ Power supply: AC 220 V \pm 10%, 50Hz
- ★ Thermometer: 0 ~ 400 °C
- ★ Ambient temperature: \leq 35 °C
- ★ Relative humidity: \leq 85%
- ★ Total power consumption: less than 650 W

★ Ordering Information:

- BGD 235----- Open Cup Flash Point Tester
- BGD 238----- Closed Cup Flash Point Tester