

# Printing Proofer

High quality proofs using gravure, gravure-offset or flexo inks are produced instantly using the BGD Printing Proofer. Featuring electronically engraved printing plates and variable printing speeds of up to 45m/min, this is an essential tool for all those involved in the manufacturing or use of liquid inks. Ideal for R & D and computer colour matching data, quality control and presentation samples. And these Printing Proofers are very easy to clean and all parts are solvent resistant.

- ◆ Almost any flexible substrate can be printed or laminated
- ◆ Excellent printability ensured by micrometer control
- ◆ Multiple ink samples may be printed simultaneously for comparison purposes
- ◆ High printing speeds enable use of inks nearing press viscosity



Ordering Information → Parameter ↓	BGD 622/1	BGD 622/2
Applicable for	Gravure Ink	Flexographic Ink
Motor Power	120W	120W
Impression roller	φ 76 × 130mm	φ 50 × 130mm
Roller Hardness	Shore A 55	Shore A 50
Printing Speed	Four sorts : 20-30-40-45m/min	Four sorts : 20-30-40-45m/min
Printing Length	150mm	150mm
Power	220V 50Hz	220V 50Hz
Overall Size	500 × 425 × 350mm	500 × 425 × 350mm
	26kg	27kg